

# Mechanical Flat Scales

Available for hire or purchase for your clinical trial.



## Overview

These mechanical flat scales provide precision, reliability and incorporate a good design.

They feature a large, clear and well defined dial making it easy to read the measurement.

They have a non-slip platform making it easy for the patient to step on and off.

## Operating Features

- Dial:
  - large area
  - well-defined
- Platform:
  - large area
  - low level
  - non-slip
  - scratch resistant
- Weight display immediately after stepping on the scale
- Reliable corrosion protection
- Galvanised frame, load lever, counterweight and all sheet metal parts

## Specification

- Capacity: 150 kg
- Graduation: 1 kg
- CE approval class: ( IIII )
- Weight: 4 kg
- Platform Dimensions: 280w x 250d mm
- Dimensions: 105h x 300w x 460d mm



**Please contact us for a quotation for the hire or purchase of this product**

Product description: - Mechanical Flat Scales

Product code - WD7PAT002

Product specifications can be subject to change without prior notice. Version: V2 10/19

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company