

## ► Mechanical Flat Scales

*available for hire or purchase for your clinical trial.*

### Overview

These mechanical flat scales with 1 lbs graduation provide precision, reliability and incorporate a good design.

They feature a large, clear and well defined dial making it easy to read the measurement.

They have a non-slip platform making it easy for the patient to step on and off.

### Specification

- Capacity: 320lbs/150 kg
- Graduation: 1 lbs, 500g
- Functions: reset-to-zero function
- Weight: 7.7 lbs / 3.5kg
- Dimensions (WxHxD): 12 x 4.7 x 18.5" / 305 x 120 x 470 mm

### Operating Features

- Dial:
  - large area
  - well-defined
- Platform:
  - large area
  - low level
  - non-slip
  - scratch resistant
- Weight display immediately after stepping on the scale
- Rugged metal housing withstands longterm use
- Galvanised frame, load lever, counterweight and all sheet metal parts



*Please contact us for a quotation for the hire or purchase of this product.*

**Product description** - Mechanical Flat Scales  
**Product code** - WD7PAT089



#### UK Office

T: +44(0)8456 777001

F: +44(0)8456 777002

E: [enquiries@woodleyequipment.com](mailto:enquiries@woodleyequipment.com)

#### USA Office

T: 1 800 471 9200 / 1 718 606 0516

F: 1 718 606 0955

[www.woodleyequipment.com](http://www.woodleyequipment.com)

**WOODLEY**  
EQUIPMENT COMPANY LTD.