

Mechanical Measuring Rod For Children & Adults

Available for hire or purchase for your clinical trial.



Overview

This mechanical measuring rod is suitable for both adults and children and can be adjusted to suit. Ideal for use in many different clinical settings.

The measurements are displayed in a window on the side of the rod and the graduation index can be held in place by the locking screw on the headpiece.

The scale is printed in centimetres on one side of the tape and centimetres and inches on the other.

Operating Features

- Simple assembly
- Easy to read
- Cost effective with reduction to the essentials
- Variable scale of 3.5 to 230 cm that is inserted into a 138 cm frame.
- Two different measurement ranges
- Wall mountable
- Choice of units of measure
- Adjustable head position

Specification

- Measuring Range: 138cm / 55.3 inches, variable from 3.5 - 230 cm / 1.4 - 90 inches
- Graduation: 1 1.7 kg / 3.7 lbs
- Dimensions: 1500h x 120w x 217d mm / 59.1h x 4.7w x 8.5d inches



Please contact us for a quotation for the hire or purchase of this product

Product description: Mechanical Measuring Rod for Children and Adults

Product code - WD7PAT099

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company