

Microman E Pipette 100 – 1000ul



Available for hire or purchase for your clinical trial.

Overview

Microman E pipette is the ideal solution to improve accuracy and precision when pipetting hot or cold liquids and is also ideal for pipetting viscous and dense samples.

The pipette works like a syringe and capillary pistons eliminate the air cushion between the sample and the disposable piston, providing 100% protection against pipette contamination and producing accurate, reproducible, and precise results every time.

The clip fitting system allows for a capillary piston to fit on to Microman E as simply as fitting a standard tip to a standard pipette.

Operating Feature

- The volume securing button allows for easy volume adjustment and minimizes the risk of accidentally changing the volume during pipetting cycles
- The Microman E pipette is comfortable and easy to use
- The shape of the handle, with an added finger hook, rests comfortably in your hand
- Fully visible volumeter on the front of the pipette so that the volume can be seen while pipetting
- A large push button allows you to aspirate and dispense with ease

Specification

- | | | | |
|--------------------------|----------|-------|-----------|
| • Volume (µL): | Min. 100 | 500 | Max. 1000 |
| • Systematic Error (µL): | ± 3 | ± 5 | ± 8 |
| • Random Error (µL): | ≤ 1.6 | ≤ 2.5 | ≤ 4 |



Please contact us for a quotation for the hire or purchase of this product

Product description -Microman E Pipette 100 – 1000u

Product Code: WD7MISC194

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company