




Microscope Eyepiece Adapter for Smartphones

Mounts most smartphones to most microscopes in seconds



CAPTURE & SHARE
MICROSCOPIC IMAGES



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"



WOODLEY
EQUIPMENT COMPANY LTD.



Microscope Eyepiece Adapter for Smartphones

InSight™ Microscope Adapter mounts most smartphones to most microscopes in seconds!

The smartphone camera can then capture photos and videos for real-time viewing on smartphone screen, via video conferencing, emailing or sharing online.



Capture and share microscopic images



50x Microscope Lens

InSight™ 50x Microscope Lens Adapter mounts onto any smartphone to capture high quality microscopic images and videos.

- Enhances the Smartphone camera
- No batteries required
- Ideal tool for
 - Dermatology Examination
 - Ectoparasite Identification
 - Teaching Aid
 - Textile Examination
 - Scientific & Laboratory Research Tool
 - Industrial Inspection
 - Computer Parts Inspection
 - Capture & Share Microscopic Images



SPECIFICATIONS

- 50 x Optical Magnification
- View Angle: Full Screen
- Focus Distance: 8-15mm
- Lens Diameter: 18.4mm
- Lens Housing: Optical Glass with Clip Holder
- Weight: 15.4g



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

WOODLEY
EQUIPMENT COMPANY LTD.