



InSight™

Microscope Eyepiece Adapter for Smartphones

**Mounts most smartphones to
most microscopes in seconds**



**CAPTURE & SHARE
MICROSCOPIC IMAGES.**



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

InSight™ Microscope Eyepiece Adapter for Smartphones

InSight™ Microscope Adapter mounts most smartphones to most microscopes in seconds!

The smartphone camera can then capture photos and videos for real-time viewing on smartphone screen, via video conferencing, emailing or sharing online.



Capture and share microscopic images

InSight™ 50x Microscope Lens

InSight™ 50x Microscope Lens Adapter mounts onto any smartphone to capture high quality microscopic images and videos.

- Enhances the Smartphone camera
- No batteries required
- Ideal tool for
 - Dermatology Examination
 - Ectoparasite Identification
 - Teaching Aid

SPECIFICATIONS

- 50 x Optical Magnification
- View Angle: Full Screen
- Focus Distance: 0.7mm
- Weight: 4.2g



+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**