

InSight

URINALYSIS

Urinalysis is an important tool in disease detection, as well as monitoring and screening animal health. Abnormalities can be indicative of diseases of the urinary system as well as other organ systems, including liver function, acid base status, and carbohydrate metabolism.

MS Reader & Printer



The Insight Urinalysis MS reader is a semi-automatic urine test strip reader with touch screen and printer to be used in conjunction with InSight MS Urine Test Strips. Simply dip the strip into the urine sample and place it on the strip loading plate. Results are displayed in 60 seconds.

- UPC results and ratio within 2 minutes.
- Microalbumin and Creatinine results and ratio within 2 minutes.



Print the results straight from the reader.

Follow us:  

Tel: +44 (0)1204 669033
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com



WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK