


MULTICHROME™

:: VETERINARY MICROBIOLOGY CULTURE SYSTEM ::

Vets can now perform simple microbiology tests in-house, quickly, easily and without the need for specialist knowledge



MultiChrome™ can test for:

- Outer ear infections
- Urinary tract infections
- Conjunctival infections
- Wound/abscess pathogens
- Upper respiratory infections
- Skin infections
- Mastitis
- Septicaemia
- Special MRSA specific

Colour Specific Detection and Identification of Common Pathogenic Bacteria in 24 hours

MultiChrome™ isolates Gram Positive and Gram Negative organisms and is designed to ignite individual, distinct colours within the bacteria allowing in-house detection and identification of pathogens at a fraction of the cost of an outside laboratory




S. agalactiae



S. uberis



S. aureus



Klebsiella



Citrobacter



E. coli

Six common mastitis bacteria

Tel: +44(0)1204 669033

Fax: +44(0)1204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

www.woodleyequipment.com

MULTICHROME™

:: VETERINARY MICROBIOLOGY INCUBATOR ::

Digital incubator with twelve assignable, built in timers.

MultiChrome™ Microbiology Incubator

- Small footprint (33h x 20.3w x 30.5d cm)
- Precise temperature control
- Digital read out
- Trays cradle all forms of dishes
- Clear door for observation
- 12 individual timers
- Audible / visible alerts
- Energy efficient, UL listed



The Microbiology Incubator is the perfect companion to the MultiChrome™ microbiology system with its adjustable controlled heat at 37°C (98°F) and patented, unique use of individual timers. No more guessing, just load, press the timer and walk away!



Each timer indicates “Green” during the 24 hour timing event and then turns “Red” with audible sound when ready to read

Precise 24 hour incubation that alerts you when completed

Tel: +44(0)1204 669033
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK