

# UltraAir Nebuliser

Available for hire or purchase for your clinical trial.



## Overview

This unit is designed for heavy-duty clinical use and is even capable to run continuously for up to 72 hours.

It operates completely silently and offers advanced output controls.

The backlight display and easy-to-read indicators make the use of this device remarkably simple.

## Operating Features

- Ultrasonic Nebuliser
- Quiet operation
- Can run continuously for up to 72 hours
- Backlit display with settings and timer
- Supplied with: Inhalation hose m (with a cuff, 70 cm). mouthpiece Set, 2 medication cups and power cord.

## Specification

- Power Consumption: 80VA
- Storage Air Pressure: 700-1.060 hPa
- Particle Size (Approx MMAD): 4.4µm
- Residual Volume: 10 ml
- Nebulisation Rate (by weight loss): 0-3 ml/min (adjustable)
- Aerosol % <5 µm: 65%
- Dimensions: 276 x 243 x 226 mm
- Weight: 4 kgs
- Power: 230V AC / 50hz



**Please contact us for a quotation for the hire or purchase of this product**

Product description: Nebuliser

Product Code: WD7PAT018

Product specifications can be subject to change without prior notice.

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company