

705CP-II Blood Pressure Monitor

Available for hire or purchase for your clinical trial.



Overview

The 705CP-II unit is specifically for professionals who wish to have a hardcopy of their measurement reading.

The unit is supplied with a compact printer that shows single reading, graph or list of all measurements in memory with an average of these readings.

This device also features PC connectivity and includes "Bio Information link Software". This program allows for easy tracking and downloading, enabling the patient to share the blood pressure history with their clinic or doctor. Please note a USB cable is not supplied. Available as an option.

Operating Features

- Fully automatic BPM with a printer for hard copies
- Intellisense technology (controlled inflation/deflation)
- Clear 3 Fold Display
- 28 readings in memory
- Supplied with std cuff, batteries, soft case and printer
- Ability to connect to PC with appropriate software and
- Clinically validated (International Protocol and BHS)
- Optional larger cuffs are available
- Compact and lightweight for ease of use
- Highly accurate algorithms
- Quiet operation, easy to use and easy to apply

Specification

- Measurement Range
 - Pressure: 0 to 299 mmHg
 - Pulse rate: 40 to 180 beats/min
- Accuracy
 - Pressure: ± 3 mmHg
 - Pulse rate: $\pm 5\%$ of reading
- Memory: 28 measurements
- Inflation: fuzzy-Logic controlled by electric pump
- Deflation: automatic pressure release valve
- Measurement Method: oscillometric method
- Pressure Detection: capacitive pressure sensor
- Power Source: 4 alkaline batteries 1.5V or AC/DC adapter
- Operating Temperature: 10°C to 40°C
- Operating Humidity: 15 to 90% RH
- Weight: 380g without batteries
- Dimensions: 71h x 177w x 115d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: 705CP-II with printer

Product code: WD7BPM003

Version: V2/19

Product specifications can be subject to change without prior notice

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company