

Omron 907 Blood Pressure Monitor

Available for hire or purchase for your clinical trial.



Overview

This innovative all-in-one Omron 907 blood pressure monitor is specifically designed for use in the clinical setting and is clinically validated.

It is designed for versatile use offering automatic oscillometric and manual auscultatory modes.

The Omron 907 is quiet, quick and simple, and is used by many professionals in a variety of settings. Measurements are displayed digitally.

Operating Features

- Clinically validated
- Oscillometric and auscultatory use combined in one unit
- Hide display function
- Deflation button (deflation control)
- Average reading mode
- Quiet operation
- Digital automatic upper arm model
- Intellisense™ technology
- Supplied with rechargeable batteries and adapter
- Environmentally friendly - rechargeable battery, washable cuff and reduced packaging
- Cuffs available in 3 sizes
- Mobile stand with basket (optional)

Please contact us for a quotation for the hire or purchase of this product

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Specification

- Measurement Range
 - Pressure: 0 to 299 mmHg
 - Pulse rate: 30 to 199 beats/min
- Accuracy
 - Pressure: ± 3 mmHg
 - Pulse rate: $\pm 5\%$ of reading
- Inflation: automatic with pumping
- Deflation: automatic by electromagnetic control valve
- Air Release: automatic rapid, by electromagnetic control valve
- Pressure Detection: electrostatic capacity semi-conductor pressure sensor
- Power Supply: AC adapter (230VAC, 50Hz, 20VA) or battery pack (4.8VDC, 6W)
- Electric Shock Protection: class II B Type
- Operating Temperature: 10°C to 40°C
- Operating Humidity: 30 to 85% RH
- Weight: Approximately 910g
- Dimensions: 203h x 139w x 131d mm



Product description: Omron 907 Blood Pressure Monitor

Product code: WD7BPM006

Product specifications can be subject to change without prior notice

Version: V2 10/19



woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company