

# Palm Style One Touch Sphygmomanometer

Available for hire or purchase for your clinical trial.



## Overview

The Riester® high-quality, palm-style sphygmomanometer has ergonomic design and is suitable for left or right-handed use.

The air release valve is regulated by finger-tip control. With tube connection at the top of the pressure gauge for unhindered blood pressure measurement.

## Operating Features

- Latex-free bulb for pumping up the cuff
- Suitable for right or left handed users
- Precision air release valve with fine adjustment and wear-free
- Microfilter protects air release valve and measuring system
- Plug-in tube connector at the top for easy handling
- Precise 64mm diameter aluminium scale
- Grey ABS-plastic casing
- Specially tempered copper-beryllium diaphragm, almost non-aging
- Supplied with std adult velcro cuff (24-32cm) and vinyl case
- Optional sizes and versions of cuffs available
- Also optionally available in a latex-free version

## Specification

- Maximum Error Tolerance: +/-3 mm Hg
- Diaphragm Pressured: Up to 600 mm HG
- Scale Diameter: 64 mm
- Scale: 0 to 300 mm Hg
- Weight: 200 g



**Please contact us for a quotation for the hire or purchase of this product**

Product description:  
Palm Style One Touch Sphygmomanometer  
Product code:WD7BPM021  
Version: V2 10/19



Product specifications can be subject to change without prior notice

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company