

Peak Flow Meter

Available for hire or purchase for your clinical trial.



Overview

These peak flow meters are designed specifically to accommodate new standards in performance and quality, they incorporate many subtle yet unique features making them easy to use and read.

The peak flow meters have dynamic characteristics which make them easier to read; it features an optional zoning system which prevents the zones being accidentally moved or tampered with.

Operating Features

- Designed specifically to meet EN13826 scale standard
- Movable asthma management colour zones that won't wash off
- Large, high visibility scale
- In-built internal check valve, avoiding cross infection concerns
- Meets and exceeds all international performance standards
- Unlike other PFM's the advanced design means dismantling is avoided
- Supplied with a 4-week recording chart

Specification

- Range: 60 - 900 l/min
- Scale Increment: 10 l/min
- Accuracy: $\pm 5\%$ or 10 l/min
- Repeatability : $\pm 2\%$
- Calibration: 100% individual
- Frequency Response: c100 Hz @ 350 l/min
- Dynamic State: critically damped @ 350 l/min
- Standard Compliance: ERS Statement 1997 (Eur Respir J 1997:9 Suppl.24) AS/NZS 4237:1994 NAEP Guidelines 1991
- Integral Check Value: expiration only



Please contact us for a quotation for the hire or purchase of this product

Product description: Peak Flow Meter

Product Code:

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company