

Pharmacy Refrigerator, 107 Litre

Available for hire or purchase for your clinical trial.



Overview

This 107 litre undercounter pharmacy refrigerator maintains an internal temperature between +2°C to +8°C and operates at +6°C for optimum storage of temperature controlled drugs and vaccines. The external digital temperature controller continuously monitors the minimum and maximum temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge.

The fridge is fan assisted to aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature.

Operating Features

- Open door sensor
- Lock
- Fixed feet
- Temperature display in the panel
- Temperature alarm
- 2 shelves
- 1 wire drawer
- Automatic fridge defrost
- Fan

Specification

- Colour: white
- Net capacity: 100 litres
- Gross capacity: 107 litres
- Door Type: hinged, reversible
- Door Style: solid
- Operating Temperature: + 2°C to +8°C
- Climate Class: N
- Voltage Supply: 230 V / 50 Hz
- Weight: 31.5 kg
- Exterior Dimensions: 830h x 495w x 580d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Pharmacy Refrigerator, 107 Litre

Product Code: WD7FRG009

Version: V3 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company