

# Pharmacy Refrigerator, 151 Litre

Available for hire or purchase for your clinical trial.



## Overview

This 151 litre solid door under counter pharmacy refrigerator maintains an internal temperature between +2°C and +8°C and the external digital temperature controller continuously monitors the temperature.

An audio-visual alarm will also activate if the temperature deviates outside of the set parameters. Fan assisted cooling around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. It is fully lockable and is provided with two keys preventing any unauthorised access.

## Operating Features

- Open door sensor
- Lock with 2 keys
- 2 adjustable front feet, 2 rear rollers
- Temperature display in the panel
- Temperature alarm
- Anti Microbial Handle
- 3 shelves
- 1 wire basket
- Digital, external controller
- Automatic fridge defrost
- Fan

## Specification

- Colour: White/Green
- Net capacity: 141 litres
- Gross capacity: 151 litres
- Door Type: Hinged, reversible
- Door Style: Solid
- Operating Temperature Range: + 2°C to +8°C
- Temperature Recording Accuracy: 0.1°C
- Climate Class: SN
- Voltage Supply: 230 V / 50 Hz
- Weight: 37 kg
- Exterior Dimensions: 845 x 595 x 635 mm HxWxD



**Please contact us for a quotation for the hire or purchase of this product**

Product description: Pharmacy Refrigerator, 151 Litre

Product code - WD7FRG010

Version: V4 10/19

Product specifications can be subject to change without prior notice.

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company