

# Pharmacy Refrigerator, 47 Litre



Available for hire or purchase for your clinical trial.

## Overview

This Pharmacy Refrigerator can be used to complement a larger medical fridge by allowing pharmaceutical items to be kept close to where they are needed, for example in a treatment room.

## Operating Feature

- Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller
- Lock and two keys
- Fan assisted
- Audible and visual high and low temperature alarms
- 3 shelf positions



## Please contact us for a quotation for the hire or purchase of this product

Product description: Pharmacy Refrigerator, 47 Litre

Product Code: WD7FRG372

Product specifications can be subject to change without prior notice.

## Specification

- Gross capacity: 47L
- Nett capacity: 30L
- Temperature Range: +2°C to +8°C (tolerance + or -2)
- Door Type: Solid
- Door Style: Hinged Reversible
- Number of shelves: 1
- Interior light (medical): No
- Auto or manual defrost: Automatic
- Refrigerant type: R600a (Extremely Green)
- Dimensions: 530 x 470 x 465mm (H x W x D)
- Weight: 20kg

±Gross m



### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [enquiries@woodleytrialsolutions.com](mailto:enquiries@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company