

Pharmacy Refrigerator, 82 Litre

Available for hire or purchase for your clinical trial.



Overview

This 82 litre, benchtop pharmacy refrigerator maintains an internal temperature between +2°C and +8°C and the external digital temperature controller continuously monitors the minimum and maximum temperatures.

An audio-visual alarm will also activate if the temperature deviates outside of the set parameters for more than 15 minutes. It has forced air cooling to ensure temperature stability and rapid temperature recovery after door openings. It is fully lockable and is provided with two keys preventing any unauthorised access.

Operating Features

- Open door sensor
- Lock
- Fixed feet
- Temperature display in the panel
- Temperature alarm
- 2 shelves
- 1 wire drawer
- Automatic fridge defrost
- Fan

Specification

- Colour: white
- Net capacity: 67 litres
- Gross capacity: 82 litres
- Door Type: hinged, reversible
- Door Style: solid
- Operating Temperature: + 2°C to +8°C
- Climate Class: N
- Voltage Supply: 230 V / 50 Hz
- Weight: 25 kg
- Exterior Dimensions: 660 x 502 x 540 mm HxWxD



Please contact us for a quotation for the hire or purchase of this product

Product description: -- Pharmacy Refrigerator, 82 Litre

Product Code: WD7FRG008

Version: V3 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company