


POISON TEST

ANTI-FREEZE ETHYLENE GLYCOL TESTING

Tests for the presence of Ethylene Glycol commonly found in "anti-freeze", lacquer, detergents and de-icers. These test strips report quantitative results of Ethylene Glycol in plasma. Providing critical information on exposure to this toxin for both Felines and Canines in as early as 30 minutes after ingestion.

*Feline & Canine
Sensitivity at
20mg/dl*

*Now available
from
Woodley*



Anti-Freeze Ethylene Glycol Test

- Feline & Canine sensitivity at 20 mg/dl
- Only 20 µl of Heparin plasma required
- Detects Ethylene Glycol approximately 30 minutes after ingestion
- Results in 10 minutes
- Sample window 12-14 hours post ingestion (varies by size, age and condition of animal)
- Room temperature storage

Tel: +44 (0) 1204 669033 (Option 1)
Fax: +44 (0) 1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHN STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK