

Hand Operated Sphygmomanometer

Available for hire or purchase for your clinical trial.



Overview

This Sphygmomanometer is a hand operated aneroid blood pressure monitor with analogue display that features a large, easy-to-read dial measuring 68mm.

The robust, corrosion-free precision manometer complies with the EN-1060 standards.

Operating Features

- Corrosion free precision manometer
- Clear 68 mm scale
- Supplied with latex-free bulb and standard adult cuff (24 -32cm)
- Comes in a vinyl case

Specification

- Measurement Range: 0 to 300 mm Hg
- Weight: 200 gms



Please contact us for a quotation for the hire or purchase of this product

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Product description: Hand Operated Sphygmomanometer

Product Code: WD7BPM017

Version: V2 10/19

Product specifications can be subject to change without prior notice



woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company