

Precisa® N Palm Sphygmomanometers

Available for hire or purchase for your clinical trial.



Overview

The precisa® N Palm Sphygmomanometer is ideal as a mobile, secondary or reserve device for any office practice and every hospital.

This sphygmomanometer is available in metal with an aluminium housing and simple tube - plug-in connection.

Suitable for right or left handed users.

Operating Features

- Latex bulb for pumping up the cuff
- Metal valve base
- Ergonomic spoon-shaped grip of stainless steel
- Single-tube connection
- Precision air release valve with fine adjustment and wear-free
- Microfilter protects air release valve and measuring system
- Specially tempered copper-beryllium diaphragm, with precision movement and almost non-aging
- Comprehensive optional selection of different sizes and versions of cuffs available
- Also optionally available in a latex-free version (when ordering enter LF for latex-free).

Specification

- Maximum Error Tolerance: +/-3 mm Hg
- Diaphragm Pressured: Up to 600 mm HG
- Linear, Tilted Aluminium Scale: Up to 600 mm Hg
- Weight: 200 g



Please contact us for a quotation for the hire or purchase of this product

Product description:
Precisa® N Palm Sphygmomanometer
Product Code : WD7BPM022
Version: V2 10/19



Product specifications can be subject to change without prior notice

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company