



## qLabs® Veterinary Coagulation Analyser

In-Practice point-of-care coagulation testing

### Features

- Requires 1 drop of non-anticoagulated whole blood (10µl)
- PT/APTT tested simultaneously on 1 test strip
- Room temperature test strip storage
- Fill test strip directly from syringe
- QC material available
- Citrate tube testing if required
- Option for Bluetooth docking station printer
- Low cost analyser and consumables
- 'Test count up' for immediate assessment
- Excellent correlation to central laboratory results

*Calibrated and validated  
specifically for use  
in cats and dogs*



### Simple Test Procedure

1. Insert test strip
2. Apply a drop of blood
3. Results displayed on LCD screen

*2 Year Warranty*

Tel: +44 (0) 1204 669033 (Option 1)  
Fax: +44 (0) 1204 669034  
Email: sales@woodleyequipment.com  
Web: www.woodleyequipment.com

**WOODLEY**  
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHN STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK



## qLabs® Veterinary Coagulation Analyser

In-Practice point-of-care coagulation testing

### Cartridge Types:

- PT/aPTT combo cartridge (12)

### Reference Intervals:

Canine:

- PT: 14 to 20 seconds
- aPTT: 94 to 123 seconds

Feline:

- PT: 13 to 20 seconds
- aPTT: 96 to 122 seconds

### Precision:

- In-house lot-to-lot and day-to-day precision <5.0 CV% (mean value)

### External QC Control:

- Lyophilised whole blood preparations
- Level I, Level II



### Helping diagnose conditions in relation to:

- Pre-operative checks for at risk animals
- Toxins and systemic diseases
- Hereditary bleeding disorders
- Excessive bleeding
- Disseminated Intravascular Coagulopathy (DIC)

### Simple Test Procedure:

1. Insert test strip
2. Apply a drop of blood
3. Results displayed on LCD screen

### Kit Contents:

- qLabs® Coagulation Analyser
- Power Adapter
- Lancets (10)
- Lancing Device
- User Manual
- Carry Case
- Data Manager Installation DVD
- Data Manager Operation Guide
- eCable

2 Year Warranty

Tel: +44 (0) 1204 669033 (Option 1)  
Fax: +44 (0) 1204 669034  
Email: sales@woodleyequipment.com  
Web: www.woodleyequipment.com

**WOODLEY**  
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHN STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK