

# RAPIDVET®-H

## BLOOD TYPING KITS - CANINE AND FELINE

In-practice tests for rapid determination of Canine DEA 1.1 and Feline A, B or AB blood types


These fast in-house blood typing card kits for cats or dogs are essential for transfusions and breeding decisions. They are fast, easy to use, economical and can be stored at room temperature.

### RapidVet®-H Feline

- For Feline A, B & AB
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5

*"Every dog & cat in need of a transfusion should be blood typed"*

Dr. Urs Giger, Dr.med.vet., Dip ACVIM, ECVIM



### RapidVet®-H Canine

- For Canine DEA 1.1
- Fast, easy to use, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Positive control well
- Available in kits of 5



*Use in conjunction with RapidVet Crossmatch Kits*

Tel: +44(0)1204 669033 (Option 1)  
Fax: +44(0)1204 669034  
Email: sales@woodleyequipment.com  
Web: www.woodleyequipment.com

**WOODLEY**  
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK