

seca CT8000i 2

Available for hire or purchase for your clinical trial.



Overview

The seca CT8000i-2 is the 2nd generation of our best-selling 12 lead ECG machine with an enhanced signal quality and newest interpretation algorithms. It is compact, portable and has a user-friendly design, that allows for quick, easy and reliable ECG acquisition. Just enter patient data, control signal quality, acquire ECG and print it with the integrated high-resolution thermal printer. Taking an ECG has never been easier.

Features

- Quick and easy-to-use ECG machine
- Compatible with all leading EMR systems - emis, SystemOne, Vision and others
- On-screen preview of ECG - review recorded ECG prior to printing, saves paper
- Outstanding signal quality for adult and paediatric ECG's due to high sample frequency
- Advanced clinical interpretation ETM algorithm - adult & paediatric
- Compact lightweight design: 1.9kg

Specification

- Display: 5" colour high resolution LCD display
- Memory: Storage of up to 100 ECGs
- Net Weight: 0.3 sts, 1.9kg
- ECG Channels: 12
- Size: 285 (W) x 189 (H) x 61 (D) mm
- Power Supply: Rechargeable batteries, mains adapter



Please contact us for a quotation for the hire or purchase of this product

Product Description: seca CT8000i 2

Product Code:

Version: V1 01/22

Product specifications can be subject to change without prior notice.



UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e:hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company