

# seca CT8000P-2

Available for hire or purchase for your clinical trial.



## Overview

The seca CT8000P-2 is a highly reliable and robust ECG machine with easy data access and fast transmission of ECG reports to your EMR system. It can receive and send data via Wi-Fi or LAN, preview all 12 ECG leads simultaneously on its generous high-resolution screen, has a large memory and a built-in A4 thermal printer. This enables you to process ECG records digitally and paperless or as an A4 printout.

## Features

- Robust and reliable: ideal for busy clinics
- Wi-Fi connectivity and seca link software included as standard – PDF ECGs automatically download to your PC following acquisition
- 8" colour high-resolution screen with anatomical electrode placement (hook-up adviser function)
- Compatible with all leading EMR systems – EMIS, SystemOne, Vision and others
- Brilliant signal quality – ideal for paediatric ECGs
- Capture and store 10 minutes of uninterrupted 12 lead data to document intermittent arrhythmias
- Lead reversal: Warnings and misplaced electrodes are immediately displayed before ECG acquisition
- Exercise ECG protocols for treadmills/ ergometers (optional extra)
- Print function
- Built-in thermal printer
- QTC interval - Ideal for Mental Health
- Preview screen feature
- ETM with ETM Sport Interpretation
- Wi-Fi or LAN
- Optimal signal quality
- QTC formulas = Bazett, Fredericia, Fremingham or Hodges

## Specification

- Display: 8" colour high-resolution screen with anatomical electrode placement
- Memory: Storage of up to 350 ECGs
- Net Weight: 0.7 sts, 4.5kg
- ECG Channels: 12
- Size: 384 (W) x 319 (H) x 90 (D) mm
- Power Supply: Rechargeable batteries, mains adapter



Please contact us for a quotation for the hire or purchase of this product

Product Description: seca CT8000P-2

Product Code:

Version: V1 01/2022

Product specifications can be subject to change without prior notice.



## UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

## USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e:hello@woodleytrialsolutions.com

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company