

Shandon Cytospin® Centrifuge

Available for hire or purchase for your clinical trial.



Overview

The Shandon Cytospin® is designed to protect the operator before, during and after sample preparation. Safety is a main priority for this centrifuge.

It allows processing of 12 specimens at one time. The lid of the centrifuge allows for one-handed opening and closing, and the lid remains locked at all times throughout the spin.

Variable speed: 200 - 2,000 RPM.

Operating Features

- Protected Program Memory: up to 23 routines for instant recall and protects from power loss
- Program key pad has "easy-touch" design
- Wipe clean control panel
- Three acceleration control buttons
- Control feature has a standby power saving mode
- Displays are also turned on and off by the standby key
- Safety Alarms:
 - out of balance
 - outside the speed tolerance
 - or if the lid is not properly locked
- As an extra measure of safety for the user, the unit will not spin if the lid is not locked
- Specimen safety alarm feature has been added to remind users at one-minute intervals to remove specimens, protect them from "air-drying" and ensure consistent results
- Supplied with Shandon Cytospin® Sealed Head: autoclavable, 12-chamber head seals completely

Specification

Speed Range: 200 - 2,000 RPM

- Timer: 1 - 99 minutes
- Power Requirements: 100-240 VAC, 50/60 Hz, 150 VA
- IEC 61010, CSA, UL and CE mark (IVDD compliant)
- Included Equipment: sealed head/rotor, marking pen, operator manual, service record book, standard workflow diagram
- Weight: 12 kg
- Dimensions (with lid open): 625h x 405w x 620d mm
- Dimensions: 240h x 405w x 620d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Shandon Cytospin® centrifuge

Product code - WD7CEN040

Version: V3 10/19

Product specifications can be subject to change without prior notice

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company