

# Shock Proof Palm Sphygmomanometer

Available for hire or purchase for your clinical trial.



## Overview

This palm style aneroid sphygmomanometer is shock-proof from a fall of up to 120cm.

Metal threaded tube connection at the top of the casing for fast cuff exchange and ergonomical operation. Fast, inflation due to optimized bulb design with integrated spoon.

Non-resistant opening of air release valve ensures precise control.

## Operating Features

- White face with blue needle
- Scale: 52 mm diameter
- Valve base and air duct with metal precision valve
- Two component ABS casting with integrated spoon
- Metal threaded tubing connector
- Non - resistant opening of air release valve ensures precise and comfortable operation
- Supplied with latex-free disinfectable one piece standard adult cuff (24 -32cm) and vinyl case

## Specification

- Maximum Error Tolerance: +/-3 mm Hg
- Weight: 173 g
- Other optional cuff sizes available:  
Small Adult Cuff 1 tube to fit 17 - 26 cm,  
Large Cuff 1 tube to fit 32 - 42 cm,  
Thigh Cuff 1 tube to fit 42 - 50 cm.



**Please contact us for a quotation for the hire or purchase of this product**

Product description:  
Shock Proof Palm Sphygmomanometer  
Product code: WD7BPM020  
Version: V2 10/19



Product specifications can be subject to change without prior notice

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company