

RAPIDPoint® 500

Blood Gas Systems

The RAPIDPoint® 500 Blood Gas Systems are designed to meet the challenges of point-of-care testing by providing fast, accurate and comprehensive test results.

- Fast, accurate, lab-quality test results in approximately 60 seconds
- Easy to use
- Maintenance free
- Automatic calibration and QC
- Built in printer and WiFi connectivity for easy transfer of results to practice management systems
- Adjustable touch screen
- Integrated bar-code reader
- Preset custom panels



Ideal for sites with the larger workload



SIEMENS
Healthineers

RAPIDPoint® 500 Tests

Blood Gas: pH, pCO₂, pO₂

Electrolytes: Na⁺, K⁺, Ca⁺⁺, Cl⁻

Metabolites: Glu, Lac

Co-oximetry: nBili, tHb, SO₂, O₂Hb, HHb, COHb, MetHb

Tel: +44 (0)1204 669033 (Option 1)
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

RAPIDPoint® 500

Blood Gas Systems



Easy Operation



1. Scan the sample using the integrated bar-code reader



2. Insert the sample



3. Analyse & receive results in 60 seconds

RAPIDPoint® 500 Tests

Blood Gas: pH, pCO₂, pO₂
Electrolytes: Na⁺, K⁺, Ca⁺⁺, Cl⁻
Metabolites: Glu, Lac
Co-oximetry: nBili, tHb, SO₂, O₂Hb, HHb, COHb, MetHb

Ideal for sites with the larger workload

• Dimensions

- Width 30cm (11.5")
- Depth 42cm (16")
- Height 55cm (21.5") display at highest position
- Weight 16.55kg (36.5lbs)
- Touchscreen: 21.1cm x 15.8cm (8.3" x 6.2")

• External Interfaces

- 3 x USB interfaces
- RS232 port
- 10 Base T Ethernet
- Barcode scanner



SIEMENS
Healthineers

Tel: +44 (0)1204 669033 (Option 1)
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK