

Table Top -20°C Freezer

Available for hire or purchase for your clinical trial.



Overview

This small table top -20°C freezer is ideal for sites where space is limited.

Our new lighter and quieter freezer is an ideal solution for the storage of samples and materials in a clinical setting. The units are energy efficient and cost effective providing the perfect solution to your refrigeration requirements.

A calibrated minimum/maximum thermometer is included as standard with rental orders (optional extra for purchased orders).

Operating Features

- Soft line design with integrated handle
- Thermostat control dial
- Reversible door hinge
- Adjustable wire shelf
- Manual defrost
- Adjustable feet
- Free standing
- Colour: white

Specification

- Freezer Capacity: 33 litres
- Voltage: 220 - 240V (Only)
- Frequency: 50Hz
- Sound Decibel Rating: 39 db (A)
- Freezing Capacity (kg/24 hours): 4
- Star Rating: 4 star
- Energy Efficiency Class: A+
- Weight: 16 kgs
- Dimensions: 492h x 472w x 450d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Table Top -20°C Medical Freeze

Product code - WD7FRZ001

Version: V3 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company