

Ultra Low Temperature Freezer ®



Available for hire or purchase for your clinical trial.

Overview

With an all stainless steel construction, this 35 litre capacity table top -85°C freezer incorporates a programmable controller with alarm. It is suitable for a variety of uses due to its compact design and noise-cancelling cooling aggregate.

The digital display shows the set temperature and the actual temperature. It has a visual and acoustic alarm and a 10mV/K output for connecting a temperature recorder.

Operating Features

- Large interior space
- Improved vacuum isolation
- Low power consumption
- Completely stainless steel casing
- Foamed form door with double isolation ring
- Stainless steel inside chamber for easy cleaning
- Ventilation grid on the back, left and right sides
- Control unit with foil keyboard
- Lockable door
- Optical and audio alarms
- RS 485 interface

Specification

- Temperature Range: -50°C to -85°C
- Control Precision: ± 1 K
- Effective Space: 35 litres
- Inside Material: stainless steel 4301
- Casing Material: stainless steel 4016
- Surrounding Temperature: +12°C to +30°C
- Power Requirements: 230V / 50Hz
- Weight: 81 kg
- Internal Dimensions: 280h x 425w x 300d mm
- External Dimensions: 540h x 580w x 765d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Table Top -80°C Medical Freezer

Product code - WD7442

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company