

iCup® Urine Drug Test Kit

Available for hire or purchase for your clinical trial.



Overview

The iCup® drug test is an all-in-one drug test cup that does all the testing for you. Once the donor has provided a specimen, the iCup® begins to test immediately. The results are ready in just seconds and can be revealed by peeling back the results label.

Because of the revolutionary design of the iCup® drug test, it can be photocopied. Once the results are clear, place the iCup® drug test side down on the photocopier using the custom designed copy template. The iCup® Drug Test detects 3 - 10 drugs instantly.

Operating Feature

- Features a temperature strip on each cup
- Two-part result form available
- Security seal
- Room temperature storage
- Built-in validity test
- Detects 3 - 10 drugs
- Fast, dark lines that make reading the test results easier
- Small, compact design makes it easy to store and handle
- Easy to use and read, donor and collector friendly
- 99% accurate
- Drug test lines start appearing within seconds
- Laboratory accurate, set to SAMSHA cutoff levels
- 18 month shelf life
- Xerox/photocopier feature for retention of drug test results

Specification

Drug Name	Abbr.	Cutoff
Amphetamine	AMP	1000 ng/mL
Barbiturates	BAR	300 ng/mL
Benzodiazepines	BZO	300 ng/mL
Cocaine	COC	300 ng/mL
Ecstasy	MDMA	500 ng/mL
Marijuana	THC	50 ng/mL
Methadone	MTD	300 ng/mL
Methamphetamine	mAMP	1000 ng/mL
Morphine	MOR	300 ng/mL
Opiates	OPI	2000 ng/mL
Oxycodone	OXY	100 ng/mL
Phencyclidine	PCP	25 ng/mL
Propoxyphene	PPX	300 ng/mL
Tricyclic Antidepressants	TCA	1000 ng/mL



Please contact us for a quotation for the hire or purchase of this product

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Product description - iCup® Urine Drug test Kit

Product specifications can be subject to change without prior notice. Version: V2 10/19

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company