

# Tube Roller

Available for hire or purchase for your clinical trial.



## Overview

These tube rollers are available as analogue or digital control. Providing a gentle rocking and rolling action for complete mixing. Available in six rollers design has a small space saving footprint, while the nine roller unit has a higher capacity.

The digital models offer variable speed from 5 to 60 rpm. The analogue models are controlled via an on/off switch and will run continuously at the fixed speed of 33rpm.

Rollers are robustly constructed and are easy to clean, with plastic rollers and a drip tray to collect spills.

## Operating Feature

- Available as either digital control or analogue control
- Ideal with all types of glass and plastic tube including those with soft rubber stoppers up to 100 mm diameter
- Suitable for mixing blood samples and viscous substances
- 6 rollers and 9 rollers available
- Variable speed and fixed speed option
- Supplied with a drip tray
- Can be used in cold rooms and in incubators

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

## Specification

- Roller Sizes: 340 x 30 mm all models
- Operating Temperature: 4°C to 60°C
- Full Specification: 15°C to 25°C
- Speed: Fixed 33 rpm (analogue models)  
- Variable 5 -60 rpm (digital models)
- Relative Humidity: 80% max.
- Amplitude: 16 mm (all models)
- Power Supply: 230V, 50Hz, 50W
- Dimensions (6 rollers): 110h x 565w x 240d mm
- Dimensions (9 rollers): 110h x 565w x 360d mm
- Weight: (6 rollers): 5.1 Kg, (9 rollers): 6.9 Kg



**Please contact us for a quotation for the hire or purchase of this product**

Product description - Tube Roller

Product code - Analogue & Digital 6 Roller and 9 Roller

Product specifications can be subject to change without prior notice.

Version: V3 10/19

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company