

Tympanic Ear Thermometer

Available for hire or purchase for your clinical trial.



Overview

Compact, elegant and versatile; this thermometer represents the highest proven accuracy for in-ear temperature measurement.

The unique positioning mode makes sure that improper positioning in the ear does not affect the accuracy of the reading.

Supplied with hard protective case and ten probe covers.

Operating Features

- Infrared technology ensures accurate readings
- Unique probe cover detector helps ensure proper hygienic use
- LCD display
- Ergonomic design
- Only 1 second measurement time
- Memory function recalls last reading
- Temperature can be displayed in Celsius and Fahrenheit
- Beeper guided operation
- Long life lithium batteries installed for immediate use
- Includes protective case
- No external influences on the readings due to hot or cold drinks or perspiration unlike traditional thermometers
- Disposable probe covers available

Specification

- Measurement Range: 34°C to 42.2°C (93.2°F to 108°F)
- Measurement Accuracy: $\pm 0.2^{\circ}\text{C}$
- Voltage: 3V DC (lithium battery CR2032)
- Power Consumption: 0.05 W
- Operating Temperature: 10°C to 40°C
- Operating Humidity: 30 to 85% RH
- Weight: 50g (including the battery)
- Dimensions: 93h x 46w x 57d mm



Please contact us for a quotation for the hire or purchase of this product

Product description - Tympanic Ear Thermometer

Product code - WD7PAT013

Product specifications can be subject to change without prior notice. Version: V2 10/19

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company