

-Basic Counter Height Freezer

2.8 Cu.Ft. / 79 L

Available for hire or purchase for your clinical trial.



Overview

This 20" wide counter height freezer is suitable for general purpose use. Capable of -20 °C freezing, with a lock for added security and a white exterior it is ideal for use in any setting.

Hospital grade cord with 'green dot' plug for added safety.

Operating Features

- Slim 20" width - 2.8 capacity inside a slim footprint
- Fully finished white cabinet allows the unit to be used freestanding
- Factory installed lock top mounted for convenient security
- Reversible door for added flexibility
- Manual defrost
- 4 clear pull out baskets
- Adjustable thermostat
- High temperature alarm
- Flat door liner
- 2 Level legs
- US Electrical Safety - ETL
- Canadian Electrical Safety - ETL-C

Specification

- Capacity: 2.8 Cu. Ft. / 79 Litres
- Temperature: -20 to -15° C / -5 to +5 °F
- Amps: 1.6
- Input Voltage: 115 V AC/60 Hz
- Weight: 57 lbs
- Dimensions External:
33.5 x 19.5 x 21.38 inches H x W x D
- Dimensions Internal:
26.36 x 14.0 x 16.75 inches H x W x D



Please contact us for a quotation for the hire or purchase of this product

Product description: Basic Counter Height 2.8 cu. ft

Product Code: WD7FRZ087

Version: V2 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company