

Under Counter Medical Laboratory Refrigerator 2.4 Cu.Ft./68 Litres

Available for hire or purchase for your clinical trial.



Overview

With 2.4 cu.ft/68 litre capacity this fan forced cooling refrigerator is ideal for use in medical, scientific pharmaceutical and laboratory facilities.

The door includes a front lock and high temperature warning alarm that also displays the current and high/low temperatures. This refrigerator includes hospital grad cord with a green dot plug.

Operating Features

- Fully furnished white cabinet
- Front mounted lock for a secure interior
- Automatic defrost
- Factory reversible door
- 2 adjustable wire shelves
- Easy to clean interior
- Thin line design
- Digital thermostat for precise temperature management
- 100% CFC Free
- US electrical safety - ETL
- Canadian electrical safety - ETL
- 2 level legs

Specification

- Capacity: 2.4 Cu. Ft. / 68 Litres
- Temperature Range: 0 °C to 10 °C / 35-46 °F
- Amps: 1.3
- Input Voltage: 115 V AC/60 Hz
- Weight: 36lbs
- Dimensions External:
26.0 x 18.5 x 17.75 inches HxWxD
- Dimensions Internal:
21.0 x 15.5 x 13.0 inches HxWxD



Please contact us for a quotation for the hire or purchase of this product

Product description: Under Counter Pharmacy Refrigerator
2.4 Cu. Ft.

Product Code: WD7FRG076

Version: V2 10/19

Product specifications can be subject to

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company