

Under Counter Pharmacy Refrigerator 4.1 Cu.Ft./116 Litres

Available for hire or purchase for your clinical trial.



Overview

With 4.1 cu.ft/116 litre capacity this refrigerator is ideal for use in medical, scientific, pharmaceutical and laboratory facilities.

This 20" wide refrigerator is designed for built in use under counters. It uses automatic defrost operation and is complete with an internal fan with gel packs for even cold air distribution throughout the refrigerator.

A high temperature alarm warns of temperature fluctuations and includes external temperature readout in the door.

Operating Features

- Flexible design allows built in freestanding
- Digital thermostats allows more precision in temperature control
- Factory installed lock provides security
- Automatic defrost
- Internal fan with gel packs
- Thin line design
- Digital thermostat for precise temperature management
- 3 wire adjustable shelves
- Internal lights
- 100% CFC Free
- US electrical safety - ETL
- Canadian electrical safety - ETL
- 2 level legs

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

Specification

- Capacity: 4.1 Cu. Ft. / 116 Litres
- Temperature Range: 0 °C to 10 °C / 35-46 °F
- Amps: 1.4
- Input Voltage: 115 V AC/60 Hz
- Weight: 64lbs
- Dimensions External:
33.5 x 19.25 x 21.38 inches HxWxD
- Dimensions Internal:
23.0 x 28.5 x 15.75 inches HxWxD



Please contact us for a quotation for the hire or purchase of this product

Product description: Under Counter Pharmacy Refrigerator
4.1 cu. ft.

Product Code: WD7FRG073

Version: V2 10/19

Product specifications can be subject to

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company