

Under Counter Pharmacy Refrigerator 5.5 Cu.Ft./155 Litres

Available for hire or purchase for your clinical trial.



Overview

This freestanding refrigerator has advanced features for use in medical, scientific, pharmaceutical and laboratory facilities.

The refrigerator has an automatic defrost operation complete with an internal fan with gel packs for even cold air distribution.

The digital thermostat offers improved temperature management. The door includes a front lock and high temperature warning alarm displaying current and high/low temperatures.

Operating Features

- Factory installed lock
- Digital thermostat
- Automatic defrost operation
- Internal fan with gel packs
- 3 wire adjustable shelves
- High temperature alarm includes external temperature readout
- Hospital grade cord with green dot plug
- 100% CFC Free
- Adjustable height
- 4 level legs
- US electrical safety - UL
- Canadian electrical safety - ULC

Specification

- Capacity: 5.5 Cu. Ft. / 155 Litres
- Temperature Range: 0 °C to 10 °C / 35-46 °F
- Amps: 1.3
- Input Voltage: 115 V AC/60 Hz
- Weight: 104lbs
- Dimensions External:
33.5 x 23.63 x 23.5 inches HxWxD
- Dimensions Internal:
27.75 x 20.0 x 17.5 inches HxWxD



Please contact us for a quotation for the hire or purchase of this product

Product description: Under Counter Pharmacy Refrigerator
.5.5 cu.ft. / 155 Litres

Product Code: WD7FRG074

Version: V2 10/19

Product specifications can be subject to

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company