

Undercounter -20°C Freezer

Available for hire or purchase for your clinical trial.



Overview

This undercounter -20°C freezer is ideal for sites where space is available under counters/worktops.

The undercounter freezer is an ideal solution to storage of samples and materials in a clinical setting. The units are energy efficient and cost effective providing the perfect solution to your refrigeration requirements.

A calibrated minimum/maximum thermometer is included as standard with rental orders (optional extra for purchased orders).

Operating Features

- Soft line design with recessed handle
- Adjustable thermostat control dial
- Reversible door hinge
- 3 freezer drawers
- Manual defrost
- Adjustable feet
- Free standing

Specification

- Freezer Capacity: 70 litres
- Voltage: 220 - 240V (Only)
- Frequency: 50Hz
- Number Of Compressors: 1
- Energy Efficiency Rating: A+
- Freezer Rating: 4*
- Sound Decibel Rating: 40 db (A)
- Weight: 34 kgs
- Dimensions: 845h x 501w x 540d mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Undercounter -20°C Medical Freezer

Product code - WD7FR7002

Version: V3 10/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company