

Undercounter -80°C Medical Freezer

Available for hire or purchase for your clinical trial.



Overview

This Ultra Low freezer can be used as either a benchtop or undercounter model. The compact space-saving design and low noise level makes it a perfect storage solution for your laboratory. This versatile, low energy consumption freezer has approximately 48 hours battery backup for acoustic alarm and temperature display in case of power cut.

Operating Features

- Single compressor technology
- VIP panels to secure lowest energy consumption
- Low energy consumption
- Low noise level
- Low heat dissipation
- Key lock on door
- Porthole for external probe
- One shelf (adjustable)
- 100% HCFC/CFC free

Specification

- Temperature Range: -40°C to -80°C
- Net Capacity: 35 litres
- Insulation: 100mm
- Casting Material: Stainless Steel 4016
- Maximum Ambient Temperature: 25°C
- Power Requirements: 230V / 50Hz (only)
- Noise Level: >45DBA
- Weight: 60kg
- Refrigerant: EP88
- Internal Dimensions: WxDxH 336x356x295 mm
- External Dimensions: WxDxH 495x660x680 mm



Please contact us for a quotation for the hire or purchase of this product

Product description: Undercounter -80°C Medical Freezer

Product Code: WD7ULT001

Version: V2 11/19

Product specifications can be subject to change without prior notice.

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company