

# Undercounter Refrigerator

Available for hire or purchase for your clinical trial.



## Overview

This undercounter refrigerator is ideal for sites where space is available under counters/worktops.

The undercounter refrigerator is an ideal solution to storage of samples and materials in a clinical setting. The units are energy efficient and cost effective providing the perfect solution to your refrigeration requirements.

A calibrated minimum/maximum thermometer is included as standard with rental orders (optional extra for purchased orders).

## Operating Features

- Soft line design with recessed handle
- Adjustable thermostat
- Internal light
- Reversible door hinge
- 2 in-door storage shelves
- 3 adjustable shelves
- Automatic defrost function
- Adjustable feet
- Free standing

## Specification

- Fridge Gross Capacity: 115 litres
- Fridge Net Capacity: 112 litres
- Voltage: 220 - 240V (Only)
- Frequency: 50Hz
- Energy Efficiency Class: A
- Sound Decibel Rating: 40 db (A)
- Weight: 29 kgs
- Dimensions: 845 x 501 x 540 mm H x W x D



**Please contact us for a quotation for the hire or purchase of this product**

Product description: Undercounter Medical Refrigerator

Product code - WD7FRG002

Version: V2 10/19

Product specifications can be subject to change without prior notice.

### UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

### USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: [hello@woodleytrialsolutions.com](mailto:hello@woodleytrialsolutions.com)

[www.woodleytrialsolutions.com](http://www.woodleytrialsolutions.com)

**woodley**  
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company