

# VISION V4002 MICROSCOPE

## HIGHLY RELIABLE AND COST-EFFECTIVE



The Vision V4002 binocular microscope is designed to fulfil all general microscopy requirements. It is a reliable, high-quality and very cost effective microscopy solution. The optical and mechanical parts of the microscope are manufactured to an exceptionally high standard, making this microscope very high value for money. It comes complete with an internal rechargeable battery and additionally a plano-concave mirror for use in the field.

### Features:

- Wide-field eyepieces - 10x/18mm
- Interpupillary adjustment 48 to 75mm
- Dioptic adjustment on left eyepiece
- 30° inclined head rotates 360°
- DIN achromatic objectives (4x, 10x, 40x, 100x)
- Variable light intensity
- LED illuminator 50,000 hours use
- Internal NiMH rechargeable battery (Up to 15 hours use from a single charge)
- N.A. Abbe 1.2 condenser
- Coaxial coarse and fine focusing
- Double layer mechanical stage - 125x116mm
- Wide range voltage for international markets – 100 to 240V
- Plano-concave mirror
- CE certified
- Product Code: WD6207



Visit [www.woodleylabdiagnostics.com](http://www.woodleylabdiagnostics.com) for full product range