



Vision V5000 Microscope

Innovative LED Powered Microscope



FIRST CHOICE

VETERINARY MICROSCOPE.

2 YEAR WARRANTY.

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.



Vision V5000 Microscope

Innovative LED Powered Microscope

The Vision V5000 Binocular Microscope has an innovative LED powered illumination system that results in a white light instead of the yellow light produced by filament lamps. It is robust and consists of high quality components which assure a high image contrast and excellent optical resolution. With a modern and ergonomic design, it is easy to handle and transport.



FEATURES

- 50,000 hours LED light source life time (instead of 1,500 from a halogen lamp)
- 30° inclined head rotates 360°
- Wide-field eyepieces - 10x/18mm
- Adjustable interpupillary distance
- Dioptic compensation
- Revolving nose piece, quadruple, reversed
- DIN achromatic objectives 4x, 10x, 40x, 100x
- Double layer mechanical stage, 125 x 115 mm with X-Y 70 x 30 mm moving range
- Coaxial coarse and fine focusing
- LED illuminator with adjustable intensity
- N.A. 1.25 Abbe condenser with adjustable height and iris diaphragm
- Universal voltage for international markets - 100 - 240V
- **Low cost UK annual service option**
- **UK lens cleaning service available free of charge**



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**