

WOODLEY EQUIPMENT COMPANY

PRODUCT GUIDE

Issue 2.5



*InSight*TM Analysis



*InSight*TM

Lab Manager



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

WELCOME TO OUR PRODUCT GUIDE

Woodley Equipment Company was **established in 1989** and prides itself on being a market leading supplier of specialist **laboratory equipment, diagnostic and critical care** products to the **veterinary industry worldwide**.

A dedicated on-the road Account Manager team provides **in-practice support**, training and on-going advice. This team is supported at our headquarters with a technical support and service department staffed with **Biomedical Scientists** and **Service Engineers** who deliver proven high quality after-sales care.

As well as a **UK nationwide sales force**, Woodley has an international distribution network offering a rapidly growing portfolio of **high quality products** with worldwide distribution rights. Woodley also has a proven track record and outstanding **customer loyalty** in the form of repeat business and client referrals.

Woodley seeks to supply veterinarians and distributors with the equipment and services to best support Practices, thus enabling them to respond to the needs of their clients, while operating a more **efficient and profitable business**.

Quantum Vet Diagnostics, a division of Woodley Equipment, is the **UK's number 1** provider of veterinary wet chemistry solutions.



"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"



epoc Critical Care/Blood Gas Analyser

The UK's leading handheld Critical Care/Blood Gas Analyser!

Test Card Contains All Parameters

Measured Analytes:

pH, pCO₂, pO₂, Na⁺, K⁺, Cl⁻, Ca⁺⁺, Glu, Lac, Hct, BUN, Creatinine, TCO₂

Calculated Values:

HCO₃⁻, BE (ecf), BE (b), sO₂, Hgb, AGapK

Features

- Accurate blood gas, electrolyte and critical care biochemistry results in 3 minutes
- Room temperature card storage
- 92µl sample, fresh whole blood (arterial, venous or capillary)
- Veterinary species specific software
- Portable handheld system with rechargeable battery & Bluetooth printer
- Ready to use anytime, anywhere: in-practice, theatre, patient side or equine yard
- Wi-Fi and Bluetooth enabled - easy transfer of data



Approved Partner

SIEMENS
Healthineers

epoc + Hemocue WBC = Rapid Pathological Screening 24 hours a day

RAPIDPoint 500e Blood Gas Systems

Ideal for sites with a larger workload

RAPIDPoint 500e Tests

Blood Gas: pH, pO₂, pCO₂

Electrolytes: Na⁺, K⁺, Ca⁺⁺, Cl⁻

Metabolites: Glu, Lac

Co-oximetry: nBili, tHb, SO₂, O₂Hb, HHb, COHb, MetHb

Features & Benefits

- Accurate, lab-quality test results in 60 seconds
- Easy to use, reliable, maintenance-free operation with cartridge based system, automatic calibration and QC
- Robust accuracy from sample to sample by leveraging the Siemens Healthineers-designed Integri-sense™ Technology

Approved Partner

SIEMENS
Healthineers



InSight™
qLabs® Veterinary
Coagulation Analyser

In-practice Point of Care coagulation testing

- Requires 1 drop of non-anticoagulated whole blood (10µl)
- PT/aPTT tested simultaneously on 1 test strip
- Room temperature test strip storage
- Fill test strip directly from syringe
- QC material available
- No need to transfer sample into citrate tube
- Low cost analyser and consumables
- 'Test count up' for immediate assessment
- 2 year warranty



HemoCue® WBC Analyser

Simple to use & results in 2 minutes

HemoCue WBC provides a total white cell count. The system offers an innovative, quick and reliable way to assess the health of the animal.

- Easy to use
- Results in 2 minutes
- Small sample volume - 10µl
- WBC count available using whole blood and synovial fluid



Partner with InSight HCT Meter for Rapid Haematological Assessment

InSight™ Veterinary
HCT Meter

Calibrated and validated specifically for Dogs, Cats and Horses

The InSight HCT Meter is a low cost handheld meter for the measurement of Canine, Feline and Equine haematocrit and total haemoglobin in capillary and venous whole blood.

- Calibrated for Dogs, Cats and Horses
- 30 second test time
- One drop of whole blood (Venous, Capillary or EDTA)
- HCT results reported as %



i-SMART 30 VET ELECTROLYTE ANALYSER



Fast & Accurate

- Fast, accurate analysis
- Tests electrolytes **Na⁺, K⁺ & Cl⁻**
- 60µl sample volume
- Results in 35 seconds
- Species specific reference ranges

All-in-One Cartridge

- Cartridge contains electrodes, all consumables & waste
- No user maintenance
- Automatic calibration

Measured electrolytes

Na ⁺	20 - 250 mmol/L
K ⁺	0.5 - 20.0 mmol/L
Cl ⁻	20 - 250 mmol/L

*Fast, accurate & maintenance free
In-Practice electrolyte analysis*

PocketChem™ BA Blood Ammonia Measurement

PocketChem BA measures blood ammonia levels. The device enables immediate testing and delivers fast, accurate results.



- Whole blood measurement
- Ready for use 24 hours a day
- Small 20µl sample volume
- Battery powered
- Easy to maintain
- Temperature compensation
- Lightweight and compact design
- Detachable printer available
- Supports the diagnosis of hepatic cephalopathy
- Helps detect and evaluate urea cycle defects
- Monitor patients with liver disease

*World's only
Point of Care
Blood Ammonia
Analyser*

InSight™ Veterinary Blood Urea Test Strips

A quick and easy, low cost test for in-practice quantitative determination of urea in whole blood.

- One drop of venous or capillary whole blood
- Test time only 90 seconds
- Colour chart for reading test results is printed on the bottle label



StatStrip® Lactate Xpress Meter

- Accurate, quantitative results in 13 seconds
- Easy to use and portable
- Small 0.6µl sample size
- Venous or arterial whole blood
- Corrects results for HCT
- No calibration required
- Eliminates interfering substances
- Stores 400 results



StatSensor® Creatinine Xpress Meter

- Rapid assessment of kidney function
- Accurate, quantitative results in 30 seconds
- Easy to use and portable
- Small 1.2µl sample size
- Capillary, venous or arterial whole blood
- Correct results for HCT
- No calibration required
- Eliminates interfering substances
- Stores 400 results



NovaVet Blood Ketone Meter


- Measures blood beta-hydroxybutrate (BHBA), the gold standard for detecting subclinical ketosis
- Accurate, quantitative results in 10 seconds
- Easy to use and portable
- No calibration required
- Simple test procedure in the field
- Small 0.8µl sample size
- Stores 400 results



Poison Test - Anti-Freeze Ethylene Glycol Testing

Tests for the presence of Ethylene Glycol commonly found in 'anti-freeze', lacquer, detergents & de-icers. These test strips report quantitative results of Ethylene Glycol in plasma.

- Feline & Canine sensitivity at 20 mg/dL
- Only 20µl of Heparin plasma required
- Detects Ethylene Glycol approximately 30 minutes after ingestion
- Results in 10 minutes
- Sample window 12-14 hours post ingestion
(varies by size, age and condition of animal)
- Room temperature storage





Veterinary Glucose Meter

g-Pet PLUS Veterinary Glucose Meter is calibrated specifically for use in Dogs, Cats and Horses, ensuring an accurate result is reported for each species.

- Biosensor auto code technology
- Accurate results
- Easy to use
- Small sample of whole blood (7µl)
- Accurate results displayed within 5 seconds
- Readable range 1.1 - 33.3 mmol/L
- 500 result automatic memory
- HCT range 20-60%
- Lifetime warranty



Calibrated for Canine, Feline & Equine use



Veterinary Lactate Meter

L-Pet only requires a small drop of blood (<3µl). The test strips draw blood in and then measure using an advanced biosensor specific for lactate concentration in blood.

- Advanced Biosensor Test Strips
- Small sample size (<3µl)
- Results in less than 45 seconds
- Large LCD display
- Stores up to 300 results
- Easy operation
- Portable
- Specific Lactate Analyser QC Material available
- 5 year warranty



An easy to use handheld meter providing blood lactate results in just 45 seconds

InSight™ Vet Vital Signs

InSight Vet Vital Signs Handheld Monitor is ideal for use in Veterinary Critical Care, General Health Screening and Anaesthesia Monitoring.

- Small handheld monitor
- Measures NIBP + SpO₂ + TEMP + Heart Rate + Mean Arterial Pressure (MAP)
- Includes cage clip & charging dock
- Easy to use
- Accurate measurement in 0.05% Low Perfusion Index
- Audio and visual alarm
- Records up to 48 hours of measurement data



Measure Blood Pressure, SpO₂, Heart Rate, Temperature & Mean Arterial Pressure (MAP)

InSight™ Vet Capnograph

InSight Vet Capnograph Monitor is ideal for use in Veterinary Critical Care and Anaesthesia Monitoring.

- Small handheld monitor
- Use to measure SpO₂ + ETCO₂ + Heart Rate
- Includes cage clip & charging dock
- Easy to use
- Accurate measurement in 0.05% Low Perfusion Index
- Audio and visual alarm
- Records up to 48 hours of measurement data



Measures ETCO₂, SpO₂ & Heart Rate

AliveCor® Veterinary ECG/Heart Monitor

Veterinary iPhone ECG and Heart Monitor is designed to snap onto your iPhone like a case and wirelessly communicates with the AliveECG Vet app. Once partnered with the ECG monitor, it transforms your iPhone into a clinical-quality, single-lead ECG recorder. Alternatively the Universal Adapter (adhesive electrode) option can attach to most smartphones.

- Easy to use - No stress or shaving required in most animals
- Powerful - Leverage the power, display and communication capabilities of your iPhone/Smartphone
- Accurate - Veterinary testing demonstrates accuracy in rate and rhythm compared to traditional devices
- Accessible - Securely store ECG notes and patient information on the device and in the cloud

Compatible With:

- iPhone 4/4S/5/5S/6/6S/SE/7 • Universal Adapter fits most Smartphones



An ECG device as mobile as you are

Vcheck V200 POINT OF CARE ANALYSER

Simplify In-Practice Point of Care Testing with Reference Lab Quality - avoiding test referrals

Achieve true analysis with fluorescence technology

- Fully quantitative automated analyser
- Accurate diagnosis for Dogs and Cats in minutes
- Analyser provides print out and result storage



Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- CRP
- Progesterone
- **Core Antibody:** Parvovirus, Distemper & Hepatitis

Feline Tests Available:

- T4
- fPL
- SAA
- **Core Antibody:** Calicivirus, Herpes & Parvovirus



Small analyser - low cost per test

Low Cost, In-Practice, Reference Lab Quality Immunoassay Testing

InSight™ Lab Manager

Results available on your smartphone via:

A stand alone Lab Management System to integrate with your Practice Management System



InSight™ Laboratory Diagnostics

BIONOTE

Anigen™ Rapid Test Kits

Available exclusively from Woodley

Feline Tests Include:

Anigen Rapid FIV Ab/FeLV Ag

Feline Immunodeficiency virus antibody and Feline Leukemia virus antigen in Feline whole blood, plasma or serum.

Anigen Rapid FCoV Ag

For the detection of Feline Corona virus antigen in Feline faeces.

Anigen Rapid FCoV Ag/Ab

Feline Corona virus antigen/antibody in Feline faeces (Ag) and whole blood, plasma or serum (Ab).

Anigen Rapid FIV Ab

Immunodeficiency virus antibody in Feline whole blood, plasma or serum.

Anigen Rapid FCoV Ab

Feline Corona virus antibody in whole blood, plasma or serum.

Anigen Rapid FPV Ag

Feline Parvo virus antigen in Feline faeces.

Anigen Rapid FHW Ag/Ab

Feline heartworm antigen and antibody in whole blood, plasma or serum.

Anigen Rapid FeLV Ag

Feline Leukemia virus antigen in Feline whole blood, plasma or serum.

Anigen Rapid FeliD-3 Ag

Feline Parvo, Corona and Giardia antigen in Feline faeces.

Anigen Rapid Toxoplasma Ab

Toxoplasma gondii antibody in Feline whole blood, plasma or serum.



Canine Tests Include:

Anigen Rapid CPV Ag

Canine Parvo virus antigen in Canine faeces.

Anigen Rapid CCV Ag

Canine Corona virus antigen in Canine faeces.

Anigen Rapid CIRDC-3 Ag

Canine Infectious Respiratory Disease (CDV, CAV and CIV) in conjunctiva (CDV Ag) and nasal fluid (CAV and CIV Ag).

Anigen Rapid CPV/CCV/Giardia Ag

Canine Parvo, Corona and Giardia antigen in Canine faeces.

Anigen Rapid Rabies Ag

Rabies virus antigen.

Anigen Rapid CPV Ab

For the titration of Canine Parvo virus in whole blood, plasma or serum.

Anigen Rapid CDV Ab

For the titration of Canine Distemper virus antibody in whole blood, plasma or serum.

Anigen Rapid CDV Ag

Canine Distemper virus antigen in conjunctiva (recommended), urine, serum or plasma.

Anigen Rapid CDV Ag/CAV Ag

Canine Distemper virus and Adeno virus antigen in conjunctiva (CDV Ag) and nasal fluid (CAV Ag).

Anigen Rapid C.Brucella Ab

Canine Brucella antibody in whole blood, plasma or serum.

Anigen Rapid Leishmania Ab

Leishmania infantum antibody in whole blood, plasma or serum.

Anigen Rapid E.Canis Ab

Ehrlichia Canis antibody in whole blood, plasma or serum.

Anigen Rapid Anaplasma Ab

Anaplasma phagocytophilum/Anaplasma platys antibody in whole blood, plasma or serum.

Anigen Rapid CPV Ag/CCV Ag

Canine Parvo virus and Corona virus antigen in Canine faeces.

Anigen Rapid Giardia Ag

Giardia antigen in Canine faeces.

Anigen Rapid CaniV-4

Dirofilaria immitis antigen, Ehrlichia canis antibody, Leishmania infantum antibody, Anaplasma, phagocytophilum/ Anaplasma platys antibody in whole blood, plasma or serum.

Anigen Rapid E.canis/ Anaplasma Ab

Ehrlichia canis and Anaplasma phagocytophilum/Anaplasma platys in whole blood, plasma or serum.

Anigen Rapid Total IgE

Canine Total IgE in serum.

This list is not exhaustive - please contact us with your requirements

NOW YOU HAVE A CHOICE

The 'Woodley Lab Systems'

Haematology, Clinical Chemistry, Immunofluorescence, Coagulation & Urinalysis Analysers

System 1



InSight™ V5 (No Retics) Haematology Analyser with Laser

Now includes InSight Lab Manager



skyla VB1 Chemistry Analyser

System 2



InSight™ 5 Diff Retic Haematology Analyser with Laser

Optional:



Vcheck V200 Immunofluorescence Analyser



InSight™ qLabs® Vet Coagulation Analyser



InSight™ Urinalysis Analyser

All your laboratory needs in one 'Complete' System for one low cost payment each month

- Includes Technical Support for 5 years
- Includes all repairs and servicing for 5 years
- Includes scheduled internal QC
- Includes 'Track your delivery'
- Includes InSight Lab Manager
- Includes software updates to improve analyser functionality and test range when appropriate
- Wet Chemistry, Endocrinology and Blood Gas options available

Contact Woodley or your Local Area Account Manager to discuss

A stand alone Lab Management System to integrate with your Practice Management System

InSight™

Lab Manager



Results available on your smartphone or tablet - access Lab Results anytime, anywhere via personal Cloud Based Files



**Wet Chemistry,
Endocrinology &
Blood Gas Options**

- InSight Lab Manager integrates your Lab Equipment with the host Practice Management System
- Alternatively, InSight Lab Manager can be used as a stand alone Lab Management System
- Results are displayed alongside a Reference Range Bar Chart for easy identification of abnormal results



Contact Woodley or your Local Area Account Manager to discuss

In-Practice Urinalysis in 24 hours simplified with


InSight™ Urinalysis

A 2-stage process available 24 hours a day

1. RapidBac™ Vet

TEST FOR BACTERIURIA

- * Negative - no bacteriuria indicated
- * Positive - proceed to Flexicult Vet if antibiotic sensitivity is required



Results in 20 minutes

RapidBac™ Vet is a lateral flow immunoassay for the detection of Gram-positive and Gram-negative bacteria in Dog and Cat urine.

- Test immediately in-practice
- Results in less than 20 minutes
- Easy to use
- Room temperature storage
- 98.5% accuracy
- No special apparatus required

Proceed with Flexicult® Vet if RapidBac™ is positive and antibiotic sensitivity is required

2. Flexicult® Vet

Chromogenic In-Practice Urine Culture & Sensitivity Test



Flexicult® Vet is a 24 hour chromogenic culture test for diagnosing UTI's in Dogs and Cats. It aids the identification and quantification of bacteria and helps select the most appropriate antibiotic for effective patient management.

Results within 24 hours

Flexicult Vet contains the following antibiotics:

- Ampicillin
- Amoxicillin/clavulanate
- Cephalotin
- Enrofloxacin
- Trimethoprim/sulfamethoxazole

Agar plate contains a chromogenic substance

Causes bacterial colonies on the agar to appear with different colours depending on the bacterium species



Examples of Flexicult Vet inoculated with *E. coli* (left) and *E. faecalis* (right)

Perform in-house without the need for specialist knowledge

Incubator required - available from Woodley



InSight™ Urinalysis

Urinalysis is an important tool in disease detection, as well as monitoring and screening animal health. Abnormalities can be indicative of diseases of the urinary system as well as other organ systems, including liver function, acid base status and carbohydrate metabolism.

InSight™ Urine Test Strips

MS-11 Urine Strips Test for:

- Urobilinogen
 - Glucose
 - Bilirubin
 - Ketone
 - Specific Gravity
 - Leucocytes
 - Blood
 - pH
 - Protein
 - Creatinine
 - Nitrate
- with urine protein/creatinine ratio



MS-2 Urine Strips Test for:

- Microalbumin
 - Creatinine
- with urine microalbumin/creatinine ratio

Refractometer

- Handheld veterinary refractometer with adjustable calibration
- Provides quick and accurate indication of vital fluid levels.
- The triple scale reads: Urine Specific Gravity, Total Serum Protein and Refractive Index
- Refractometer equipped with ATC (Automatic Temperature Compensation)




InSight™ Urinalysis Reader

- Semi-automatic urine test strip reader with touch screen & printer to be used with InSight Urine Test Strips
- Simply dip the strip into the urine sample and place it on the strip loading plate
- Results displayed in 60 seconds



Urinalysis Atlas

- Sample collections
- Physical characteristics
- Examination of sediment
- Photographic identification of casts
- Contaminants in urine sediment



Use in conjunction with RapidBac Vet & Flexicult Vet for complete In-Practice Urinalysis

InSight™ Microscope Eyepiece Adapter for Smartphones

- Mounts any smartphone to any microscope in minutes!
- Smartphone camera can then capture photos & videos for real-time viewing on smartphone screen, via video conferencing, emailing or sharing online

Capture & share microscopic images

Vision V5000 Microscope

- LED illuminator
- Excellent optical resolution
- Wide-field eyepieces
- Dioptic compensation
- DIN achromatic objectives 4x, 10x, 40x, 100x
- Coaxial coarse and fine focusing



Haematology Atlas

Haematology Atlas written exclusively for Woodley Equipment by Kostas Papasouliotis DVM PhD DipRCPath MRCVS (Bristol Vet School).

Titles Available - Canine, Feline & Equine
These atlases are an excellent microscope-side reference guide offering normal cell morphology.

- Erythrocytes
- Neutrophils
- Lymphocytes
- Eosinophils
- Basophils
- Platelets
- Methodology: Microhaematocrit
- Methodology: Blood Smears



ImmunoComb® VacciCheck® In-Practice Antibody Test Kits

ImmunoComb VacciCheck Antibody Test Kit monitors serological status, reduces vaccination failure and prevents over vaccination consequences for Dogs and Cats.

VacciCheck is recommended by the WSAVA to comply with regulations to reduce unnecessary vaccinations

VACCICheck® Canine

The Canine VacciCheck Antibody Test Kit is designed to determine Dog serum antibody titre to:

- Infectious Canine Hepatitis (ICH)
- Canine Parvovirus (CPV)
- Canine Distemper Virus (CDV)

The main purpose of this kit is to provide a useful tool to assess the immune status of Dogs to the three core Canine infections.



VACCICheck® Feline

Quick & easy to use


The Feline VacciCheck Antibody Test Kit is designed to determine Cat serum IgG antibody titre to:

- Feline Panleukopenia (FPLV)
- Feline Herpes Virus (FHV)
- Feline Calici Virus (FCV)

The main purpose of this kit is to provide a useful tool to assess the immune status of Cats to the three core Feline infections.



WSAVA Guidelines state that titre testing for the 3 core Canine and Feline pathogens should be performed instead of automatic annual revaccination



Clinispin MPC Combi Veterinary Centrifuge



Includes MPC Combi Rotor for 18 x HCT/QBC tubes & 4 x 2ml tubes



The Clinispin MPC Combi is a high speed centrifuge equipped with a maintenance free drive, digital display and simple interface for silent & efficient operation. The combination rotor provides efficient separation of blood and urine samples and accurate Haematocrit determination.

- 3 pre-defined presets for easy and rapid operation at one press of a button:
 - Serum/Plasma (12,000 RPM for 2 minutes)
 - HCT/QBC® (12,000 RPM for 5 minutes)
 - Urine (1,000 RPM for 1 minute)
- Presets can be varied for each centrifugation if required
- Rotor allows for quick and easy sample loading
- Maintenance free drive
- Digital display
- Easy to use
- Quiet operation
- Displays RPM or RCF (g-force)

RapidVet® - H Blood Typing Kits - Canine & Feline

In-practice tests for rapid determination of Canine DEA 1.1 and Feline A, B or AB blood types.

These fast in-house blood typing card kits for Dogs and Cats are essential for transfusion and breeding decisions. They are fast, easy to use, economical and can be stored at room temperature.


RapidVet-H Canine

- For Canine DEA 1.1
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5



RapidVet-H Feline

- For Feline A, B & AB
- Fast, easy, economical and VITAL
- Results in 2 minutes
- Auto-agglutination screen
- Requires 3 drops of whole blood
- Available in kits of 5



"Every Dog & Cat in need of a transfusion should be blood typed" Dr. Urs Giger, Dr.med.vet., Dip ACVIM, ECVIM

RapidVet® - H Crossmatch Kits

In-practice crossmatch test kits for use in transfusion medicine on either Canine or Feline patients.

- Will detect all major Canine and Feline blood group incompatibilities
- Performed in-practice in 10 minutes
- Easy to use and results are clear
- Available in kits of 3 tests
- Minor Crossmatch Kits also available

Use alongside Canine & Feline Blood Typing Kits

Essential for good transfusion practice



Veterinary Parasep® Faecal Filters

- Determine worm egg counts in Cattle, Sheep and Horses
- Cleaner, safer method
- Separates eggs from faecal filter in safe, enclosed environment
- No possibility of cross contamination
- More accurate and consistent results
- No laborious, time consuming 'clean up' of sinks, shaker bottles and sieves



EZ-SmartTemp™ Min Max Refrigeration Thermometer

- Easy to use - simple and small design with easy to read digital display
- This calibrated thermometer attaches to the exterior of the refrigeration unit and measures the internal temperature via a probe which is sealed in a non toxic glycol solution
- Bottle insulates sensor from rapid temperature changes when the door is opened
- Triple display simultaneously shows the minimum, maximum and current temperature
- Alarm feature alerts if temperature rises above or falls below a set point - programmable in 1°C increments
- Supplied with velcro and magnetic strips for mounting
- Supplied with calibration certificate



EZ-SmartLog™ One USB Temperature Datalogger

- Measures and stores up to 16,382 temperature readings over a -35°C to +80°C range
- Direct USB connection for easy set-up, functionality and data download
- User-programmable alarm thresholds
- Status indication via red, green and orange LEDs
- Complete with Windows control software which is easy to download and install
- Supplied with replaceable lithium battery - typically allows logging for 1 year



InSight™ DB



Integrates with Practice Management Systems

A fully automated clinical chemistry analyser which uses liquid reagents and is approximately 70% cheaper to run than dry chemistry

Chemistry Analyser

- User friendly intuitive veterinary software
- Onboard refrigerated reagent/sample compartment
- Minimal user maintenance
- Fast throughput - example batch times:
9 mins 1x GHP, 13mins 2x GHP, 16mins 3x GHP
- Continuous loading, no need to wait for batch completion
- Stat function for urgent tests
- Fully automated with auto-reflex dilution of over range results
- Flexible testing, single parameter or up to 20 user defined profiles
- Small footprint W 414mm x D 565mm
- Online remote assistance directly to your analyser
- Long life ceramic syringes to ensure accuracy and precision
- Environmentally friendly semi-disposable reaction cuvettes with high efficiency wash system

Customer defined profiles and/or tests can be programmed from the following:

Albumin, Alkaline Phosphatase, ALT, AST, Amylase, Bilirubin (Total), Calcium, Cholesterol, CK, Creatinine, Glucose, GGT, Phosphate, Protein (Total) & Urea.

BHB, Bile Acids, Chloride, Copper, Fibrinogen, Fructosamine, GLDH, GSH-Px (as Selenium), Haemoglobin, IgG, LDH, Lipase, Magnesium, NEFA, SAA, T4 (Total Thyroxin), Triglycerides, Uric Acid, Urine Creatinine, Zinc, Phenobarbitone, Iron, Ferritin and TIBC

Fully supported Nationwide by our Specialist Team of Field Service Technicians

InSight™ p-LINK OR InSight™ Lab Manager

Connectivity Choices to link your Analysers to your Practice Management System



Cloud based access to Lab Reports

- **InSight p-LINK** is a small computer (Raspberry Pi) which enables sample results from in-practice analysers to be transferred to a Vet Practice network folder with results in Text, CSV or VetXML formats
- Now includes cloud based access to lab reports
- Quick and easy to set up via a wired LAN point or Wi-Fi
- **InSight Lab Manager** is an advanced, fully integrated Laboratory Management System with specialist IT support which integrates in-practice analysers with the patient record, ensuring results aren't lost and are readily available for diagnosis
- Advanced graphical comparison allows specific analyte measurements to be plotted, emailed and printed



AIA-360 Analyser

Offers the automation powers and performance of a large analyser in a benchtop format. Random access allows any combination of tests, in any order at any time and a throughput of 36 results per hour with the first test result in 20 minutes ensures a fast reliable service.

- Simple to use
- Fast, reliable results
- Simple maintenance
- Touchscreen with intuitive menus

Key Assays Available:

- Free T4
- Total T4
- Cortisol
- Progesterone
- ACTH
- Insulin
- cTnl

Fully supported Nationwide by our Specialist Team of Field Service Technicians

Integrates with Practice Management Systems

Take a 360° look at your endocrinology testing with the TOSOH AIA-360

InSight™ V5

A fully automated
multi-species
Haematology Analyser

**5 Part Differential
with Laser Technology**

InSight™ 5 Diff Retic

A true 5 part differential
laser Haematology Analyser
including **Reticulocytes**

**Uses dual laser and
impedance technology**

Haematology Analysers

- Quick and easy to use
- Small sample volume - 15µl whole blood
- 60 second analysis time
- 5-part differential, 23 parameters
- Species specific reference ranges
- Large user-friendly colour touch screen
- Low maintenance
- Fully automated
- Stores up to 50,000 results
- Dog, Cat, Rabbit, Horse & Cow Settings

**Both analysers integrate with
Practice Management Systems**

- Compact, benchtop size
- Small sample size - 20µl whole blood
- 5-part differential, 27 parameters
- Reticulocytes on every sample
- 2 minute cycle time
- Quick and easy operation
- Low maintenance
- Touch screen
- Over 50 species on board



skyla VB1 Veterinary Chemistry Analyser

For both Rotor Profiles
& Single/Dual Tests

**Highly accurate test results
in under 15 minutes**

Integrates with InSight Lab
Manager - a stand alone
Lab Management System
which can integrate to your
Practice Management System

- Fast and reliable
- Compact, benchtop size
- Performs up to 15 biochemistry tests simultaneously
- Performs up to 3 single/dual tests simultaneously
- 200µl whole blood, serum or plasma
- LCD colour touch screen & built-in printer
- Fully automatic system with walk away operation
- Built-in QC and calibration for each test
- Species specific reference ranges
- Stores up to 50,000 results



**Integrates with Practice
Management Systems**

InSight™ Lab Manager

VB1 Panel Test Portfolio:

- | | |
|------------------|-------------------|
| • Preanaesthetic | • Large Animal |
| • Renal | • Liver Plus |
| • Diagnosis | • Avian & Reptile |
| • Liver | • Diabetes |
| • Diagnosis Plus | • Electrolytes |
| • Equine | • Diagnosis II |
| • Thyroid/T4 | • Critical Care |

Parameters Available:

ALB, ALP, ALT, AMY, AST, BUN, Ca, Cl⁻, CREA, CPK, GGT, GLU, K⁺, Na⁺, PHOS, TBIL, CHOL, TP, B Acids, T4, Mg, UA, TCO₂, TRIG, LAC, FRU, LIPA, UCRE, UPRO, GLOB*, A/G*, Na⁺/K⁺, BUN/CREA*, UPC*, UREA* * Calculated

Single Test Dual Test

- | | |
|-------------------|------------------------------------|
| • CREA | • Na ⁺ & K ⁺ |
| • Cl ⁻ | • PHOS & Ca |
| • AMY | • CHOL & TRIG |
| • BUN | • GGT & CPK |
| • TBIL | • ALT & ALP |
| • LAC | • ALB & TP |
| | • AST & CPK |
| | • LIPA & AMY |
| | • UPC Kit |
| | (UPRO & UCRE) |



InSight™ Analysis

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

Delivering excellence to our customers for over
30
years



Approved
Partner

SIEMENS
Healthineers



01204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

WOODLEY
EQUIPMENT COMPANY LTD.