iTherm CH25 NXT

Cooling and Heating Incubator





RELIABLE TEMPERATURE,
WIDE RANGE OF CONTROL

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

iTherm CH25 NXT

Cooling and Heating Incubator

With the iTherm CH25 NXT, you get a steady and reliable temperature with a wide range of control from 15°C below ambient to 60°C. It comes with an internal power outlet for a small shaker or mixer and can hold various volumes up to 2L.

FEATURES



- Temperature control range from 15°C below ambient to 60°C
- Save sample from cross contamination
- Temperature accuracy of ±3.5°C (at 37°C at centre point without load)
- Rapid heating and uniform temperature throughout
- Large capacity of 45 litres, allows storage of large batches
- Brings down the temperature and eliminates the need for a compressor
- Internal power outlet for installing other devices
- Comes with a built-in calibration function of up to ±5°C

SPECIFICATIONS

Temperature Range:	2°C to 60°C
Controlled Temperature Range:	15°C below ambient to 60°C
Temperature Accuracy:	±3.5°C (at 37°C at centre point without load)
Temperature Uniformity:	±3.5°C (at 37°C at centre point without load)
Temperature Increments:	±1
Exterior Dimensions:	400 x 500 x 510 mm
Interior Dimensions:	325 x 325 x 425 mm
Capacity:	45L
Weight:	15kg
Electrical:	110 to 240 VAC

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY EQUIPMENT COMPANY LTD.